

## Certificate of Analysis

No. CA07178

Product: Eichrom's LN Resin

Lot #: LNA11267

Eichrom's LN Resin, Lot # LNA11267, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

| <u>Test Value</u>                     | <u>Specification</u> | <u>Result</u> |
|---------------------------------------|----------------------|---------------|
| Appearance                            | White beads          | Pass          |
| Distribution Ratio, 0.1 N nitric acid | >700                 | 1004          |



January 4, 2008

Jill Bryant  
Quality Control Chemist