

Certificate of Analysis


No. CA12018

Product: Eichrom's LN Resin
Lot #: LNS011412

Eichrom's LN Resin Lot # LNS011412 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	1038



January 27, 2012

Jill Bryant
Quality System Coordinator