

Certificate of Analysis

No. CA15052

Product: Eichrom's LN Resin

Lot #: LNS031515

Eichrom's LN Resin Lot # LNS031515 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2045.1



Sarah Teichma
Quality Manager

25 Mar 15

March 25, 2015