

Certificate of Analysis

No. CA08242

Product: Eichrom's MnO₂ Resin

Lot #: MNA05106


Eichrom's MnO₂ Resin Lot # MNA05106 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery	≥ 95%	97 %

cpm/mL vs. mL of Eluate

1.5 grams MnO₂ Resin, slurry packed from Deionized water, Flow rate = 10 mL/min, 22(1) °C



Jill Bryant

December 1, 2008

Jill Bryant
Quality System Coordinator