



eichrom
TECHNOLOGIES

A GCI COMPANY

Certificate of Analysis

No. CA17110

Product: Eichrom's MnO₂ Resin

Lot #: MNA062217

Eichrom's MnO₂ Resin Lot # MNA062217 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	≥ 95%	98.7 %
Ba-133 in Load solution (4 reps)	<1.0%	0.89 %

Sarah McAlister
Quality Manager

July 7, 2017