



Certificate of Analysis

No. CA14158

Product: Eichrom's MnO₂ Resin

Lot #: MNA091714

Eichrom's MnO₂ Resin Lot # MNA091714 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	$\geq 95\%$	102.0 %
Ba-133 in Load solution (4 reps)	$< 1.0\%$	0.00 %



Sarah Teichma
Quality Manager

24 Sep 14

September 24, 2014