

Certificate of Analysis


No. CA10169

Product: Eichrom's MnO_2 Resin cartridges
Lot #: MNAR10B

Eichrom's MnO_2 Resin cartridge Lot # MNAR10B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Ba-133 recovery	$\geq 95\%$	99.6 %



Jill Bryant

August 27, 2010

Jill Bryant
Quality System Coordinator