

Certificate of Analysis

No. CA11120

Product: Eichrom's MnO₂ Resin cartridges


Lot #: MNAR11A

Eichrom's MnO₂ Resin cartridge Lot # MNAR11A has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	≥ 95%	105.4 %
Ba-133 in Load solution (4 reps)	<1.0%	0.0%

Jill Bryant
Quality System Coordinator



May 16, 2011