

Certificate of Analysis

No. CA14015


Product: Eichrom's MnO₂ Resin cartridges

Lot #: MNAR14B

Eichrom's MnO₂ Resin cartridge Lot # MNAR14B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	≥ 95%	101.2 %
Ba-133 in Load solution (4 reps)	<1.0%	0.00 %



Sarah Teichma
Quality Manager

14 Feb 14

February 14, 2014