

Certificate of Analysis
No. CA 17041

Product: Eichrom's MnO2 Resin Cartridges
Lot #: MNAR17A

Eichrom's MnO2 Resin cartridge Lot # MNAR17A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, MNA060916, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	> 95%	99.1 %
Ba-133 in Load solution (4 reps)	<1.0%	0.91 %

Jill Bryant
Quality System Coordinator



March 17, 2017