

## Certificate of Analysis


No. CA17119

Product: Eichrom's MnO<sub>2</sub> Resin Cartridges  
Date of Manufacture: 20-Jul-17  
Lot #: MNAR17C

Eichrom's MnO<sub>2</sub> Resin cartridge Lot MNAR17C has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, MNA060916, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Ba-133 recovery (4 reps)	≥ 95%	99.1 %
Ba-133 in Load solution (4 reps)	<1.0%	0.91 %

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

July 20, 2017