

Certificate of Analysis

No. CA07084


Product: Eichrom's MnO₂ Resin Cartridges

Lot #: MNAR7B

Eichrom's MnO₂ Resin cartridge Lot # MNAR7B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
% Solids	45-60%	48.1%
Dw Ba (1 min)	>50	55
Dw Ba (3 min)	>90	99
Dw Ba (5 min)	>150	260
Dw Ba (10 min)	>250	900
Dw Ba (20 min)	>400	1900
Dw Ba (60 min)	>1000	6600



June 25, 2007

Jill Bryant
Quality Control Chemist