



# Certificate of Analysis

No. CA 14079

Product: Eichrom's Nickel Resin columns  
Lot #: NI14B

Eichrom's Nickel Resin column Lot # NI14B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	98.6 %
Ave. % Ni in Load/Rinse	<1%	0.0 %
Ave. % Ni in Strip	>99%	99.8 %

Jill Bryant  
Quality System Coordinator

May 12, 2014