

Certificate of Analysis

No. CA08011

Product: Eichrom's Nickel Resin
Lot #: NIA01038

Eichrom's Nickel Resin Lot # NIA01038 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	99.9 %
Ave. % Ni in Load/Rinse	<1%	0.1%
Ave. % Ni in Strip	>99%	99.9%



January 23, 2008

Jill Bryant
Quality Control Chemist