

## Certificate of Analysis

No. CA06043


Product: Eichrom's Nickel Resin

Lot #: NIA03296

Eichrom's Nickel Resin Lot # NIA03296 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101 %
Ave. % Ni in Load/Rinse	<1%	0
Ave. % Ni in Strip	>99%	100%



Jill Bryant  
Quality Control Chemist

Prepared: April 14, 2006