

Certificate of Analysis

No. CA11115

Product: Eichrom's Nickel Resin
Lot #: NIA042311

Eichrom's Nickel Resin Lot # NIA042311 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	102 %
Ave. % Ni in Load/Rinse	<1%	0 %
Ave. % Ni in Strip	>99%	100 %



Jill Bryant
Quality System Coordinator

May 13, 2011