

## Certificate of Analysis

No. CA07099

Product: Eichrom's Nickel Resin  
Lot #: NIA07097

Eichrom's Nickel Resin Lot # NIA07097 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101 %
Ave. % Ni in Load/Rinse	<1%	0
Ave. % Ni in Strip	>99%	100%



Jill Bryant  
Quality Control Chemist

Prepared: July 19, 2007