



Certificate of Analysis

No. CA 14124

Product: Eichrom's Nickel Resin

Lot #: NIA071714

Eichrom's Nickel Resin Lot # NIA071714 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	97.7 %
Ave. % Ni in Load/Rinse	<1%	0.4 %
Ave. % Ni in Strip	>99%	99.6 %


Sarah Tejchma
Quality Manager

28 July 2014
July 28, 2014