

Certificate of Analysis

No. CA09090


Product: Eichrom's Nickel Resin cartridges

Lot #: NIAR09A

Eichrom's Nickel Resin cartridge Lot # NIAR09A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	103 %
Ave. % Ni in Load/Rinse	<1%	0.1%
Ave. % Ni in Strip	>99%	99.9%



July 13, 2009

Jill Bryant
Quality System Coordinator