

Certificate of Analysis


No. CA09095

Product: Eichrom's Nickel Resin cartridges
Lot #: NIAR09B

Eichrom's Nickel Resin cartridge Lot # NIAR09B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101 %
Ave. % Ni in Load/Rinse	<1%	0.1%
Ave. % Ni in Strip	>99%	99.9%



August 3, 2009

Jill Bryant
Quality System Coordinator