

Certificate of Analysis

No. CA10226

Product: Eichrom's Nickel Resin cartridges
Lot #: NIAR10A

Eichrom's Nickel Resin cartridge Lot # NIAR10A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101 %
Ave. % Ni in Load/Rinse	<1%	0.1 %
Ave. % Ni in Strip	>99%	99.9 %



October 22, 2010

Jill Bryant
Quality System Coordinator