

Certificate of Analysis

No. CA 13146


Product: Eichrom's Nickel Resin cartridges
Lot #: NIAR13C

Eichrom's Nickel Resin cartridge Lot # NIAR13C has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	97.1 %
Ave. % Ni in Load/Rinse	<1%	0.3 %
Ave. % Ni in Strip	>99%	99.7 %

Jill Bryant
Quality System Coordinator



September 13, 2013