

Certificate of Analysis


No. CA07038

Product: Eichrom's Nickel Resin Cartridges
Lot #: NIAR7A

Eichrom's Nickel Resin cartridges Lot # NIAR7A has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101 %
Ave. % Ni in Load/Rinse	<1%	0
Ave. % Ni in Strip	>99%	100%



Jill Bryant
Quality Control Chemist

Prepared: March 13, 2007