



Certificate of Analysis

No. CA 14225

Product: Eichrom's Pb Resin columns

Lot #: PB14A

Eichrom's Pb Resin column Lot # PB14A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.00 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	90.4 %
Chemical Yield	80-120 %	115.6 %

Jill Bryant
Quality System Coordinator

December 1, 2014