

Certificate of Analysis

No. CA 08246

Product: Eichrom's PB Resin


Lot #: PBA11118

Eichrom's PB Resin Lot # PBA11118, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.01 %	0.00 %
% Lead in Second Rinse	< 1.0 %	0.00 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	88.0 %
Chemical Yield	80-120 %	91.5 %

Jill Bryant
Quality Control Chemist



December 2, 2008