

## Certificate of Analysis


No. CA10080

Product: Eichrom's PB Resin cartridges  
Lot #: PBSR10A

Eichrom's PB Resin cartridge Lot # PBSR10A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Gravimetric Yield	80-120 %	82.7 %
Corrected Pb Recovery	80-120 %	112 %



April 30, 2010

Jill Bryant  
Quality System Coordinator