

## Certificate of Analysis


No. CA 10254

Product: Eichrom's PB Resin cartridges  
Lot #: PBSR10B

Eichrom's PB Resin cartridge Lot # PBSR10B, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.00 %
% Lead in Second Rinse	< 1.0 %	0.00 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	89.9 %
Chemical Yield	80-120 %	119 %



Jill Bryant  
Quality System Coordinator

December 14, 2010