

## Certificate of Analysis

No. CA 11121

Product: Eichrom's PB Resin cartridges

Lot #: PBSR11B

Eichrom's PB Resin cartridge Lot # PBSR11B, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.00 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	86.6 %
Chemical Yield	80-120 %	100.1 %

  
Jill Bryant  
Quality System Coordinator

May 17, 2011