

Certificate of Analysis

No. CA 08080


Product: Eichrom's PB Resin Cartridges
Lot #: PBSR8A

Eichrom's PB Resin cartridges, Lot # PBSR8A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.01 %	0.00 %
% Lead in Second Rinse	< 1.0 %	0.00 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	99.8 %
Chemical Yield	80-120 %	93.2 %

Jill Bryant
Quality Control Chemist



April 21, 2008