

Certificate of Analysis


No. CA10145

Product: Eichrom's Pre-Filter Resin cartridges
Lot #: PFLR10A

Eichrom's Pre-Filter Resin cartridges, Lot PFLR10A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Flow Rate	> 0.80 ml/min.	Pass



July 28, 2010

Jill Bryant
Quality System Coordinator