

Certificate of Analysis

No. CA08109

Product: Eichrom's RE Resin columns


Lot #: RE8A

Eichrom's RE Resin column Lot # RE8A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 125	185

Signature



June 10, 2008

Jill Bryant

Quality System Coordinator