

## Certificate of Analysis

No. CA11097

Product: Eichrom's RE Resin

Lot #: REA041411

Eichrom's RE Resin Lot # REA041411, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 125	157

  
Sarah Tejhema  
Quality Supervisor

4-25-11

April 25, 2011