

Certificate of Analysis

No. CA12128


Product: Eichrom's RE Resin cartridges
Lot #: RESR12A

Eichrom's RE Resin cartridge Lot # RESR12A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	> 125	159.7

Jill Bryant
Quality System Coordinator



July 9, 2012