

# Certificate of Analysis

No. CA 14176


Product: Eichrom's 25mm Resolve™ Filter Funnels

Lot #: RFDF081814

Eichrom's 25mm Resolve™ Filter Funnels Lot # RFDF081814 has been manufactured according to Eichrom standard operating procedures. The filter membrane used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Am-241 FWHM avg keV (8 replicates)	<50.0	41.8
Am-241 FWHM values above 60 keV	0	0
Am-241 FWHM values between 50-60 keV	2 or fewer	0
Am-243 FWHM avg keV (8 replicates)	<50.0	40.2
Am-243 FWHM values above 60 keV	0	0
Am-243 FWHM values between 50-60 keV	2 or fewer	0

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

September 26, 2014