

## Certificate of Analysis

No. CA07008

Product: Eichrom's Sr Resin

Lot #: SRA01167

Eichrom's Sr Resin, Lot # SRA01167 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	119
Gravimetric Recovery	$\leq 105\%$	100 %
Chemical Recovery, Corrected	90-110 %	98.2 %



Prepared: January 31, 2007

Jill Bryant  
Quality Control Chemist