

Certificate of Analysis

No. CA06045


Product: Eichrom's Sr Resin

Lot #: SRA04076

Eichrom's Sr Resin, Lot # SRA04076 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	111
Gravimetric Recovery	$\leq 105 \%$	98.3 %
Chemical Recovery, Corrected	90-110 %	98.9 %



Jill Bryant

Quality Control Chemist

Prepared: April 18, 2006