

Certificate of Analysis
No. CA 15102


Product: Eichrom's Sr Resin
Lot #: SRA060515

Eichrom's Sr Resin Lot # SRA060515 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, Sr in 3.1 N nitric acid	≥ 100	126.2
% Gravimetric Recovery	≤ 105	100.8
% Chemical Recovery, Corrected	90-110	95.6
% Sr stripped in 0.05M HNO ₃		
Least acceptable mean (5 rep)	85.0	95.7
Max. standard deviation	2.5	0.7
%Ca stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.15	0.0
Max. standard deviation	0.05	0.0
%Ba stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.50	0.2
Max. standard deviation	0.50	0.0
%Y stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.0
Max. standard deviation	0.25	0.0
%Y in 8M HNO ₃ load and rinse		
Least acceptable mean (5 rep)	99.0	99.9
Max. standard deviation	2.5	0.0


Sarah Tejehma
Quality Manager


June 11, 2015