

## Certificate of Analysis

No. CA07139

Product: Eichrom's Sr Resin


Lot #: SRA08027

Eichrom's Sr Resin Lot # SRA08027 has been manufactured according to Eichrom standard operating procedures. The resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

| <u>Test Value</u>                     | <u>Specification</u> | <u>Result</u> |
|---------------------------------------|----------------------|---------------|
| Appearance                            | White beads          | Pass          |
| Distribution Ratio, 3.1 N nitric acid | $\geq 100$           | 121           |
| Gravimetric Recovery                  | $\leq 105 \%$        | 99 %          |
| Chemical Recovery, Corrected          | 90-110 %             | 97.7 %        |

Jill Bryant  
Quality Control Chemist



Prepared: September 24, 2007