

Certificate of Analysis

No. CA06088


Product: Eichrom's Sr Resin

Lot #: SRA08046

Eichrom's Sr Resin, Lot # SRA08046 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	132
Gravimetric Recovery	$\leq 105\%$	104 %
Chemical Recovery, Corrected	90-110 %	97.9 %



Jill Bryant

Quality Control Chemist

Prepared: August 14, 2006