

## Certificate of Analysis


No. CA06099

Product: Eichrom's Sr Resin  
Lot #: SRA09226

Eichrom's Sr Resin, Lot # SRA09226 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	113
Gravimetric Recovery	$\leq 105\%$	102 %
Chemical Recovery, Corrected	90-110 %	94.5 %



Jill Bryant  
Quality Control Chemist

Prepared: October 9, 2006