

Certificate of Analysis

No. CA05120


Product: Eichrom's Sr Resin

Lot #: SRA11145

Eichrom's Sr Resin, Lot # SRA11145 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	129
Gravimetric Recovery	$\leq 105\%$	101 %
Chemical Recovery, Corrected	90-110 %	94.3 %


Jill Bernstein

Quality Control Chemist

Prepared: November 29, 2005