

Certificate of Analysis

No. CA05129


Product: Eichrom's Sr Resin

Lot #: SRA12085

Eichrom's Sr Resin, Lot # SRA12085 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	138
Gravimetric Recovery	$\leq 105\%$	97.9 %
Chemical Recovery, Corrected	90-110 %	96.9 %


Jill Bernstein

Quality Control Chemist

Prepared: December 20, 2005