

Certificate of Analysis

No. CA07096

Product: Eichrom's Sr Resin

Lot #: SRF06277

Eichrom's Sr Resin, Lot # SRF06277 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	124
Gravimetric Recovery	$\leq 105\%$	103 %
Chemical Recovery, Corrected	90-110 %	94.7 %



July 10, 2007

Jill Bryant
Quality Control Chemist