

Certificate of Analysis

No. CA07004


Product: Eichrom's Sr Resin

Lot #: SRS01127

Eichrom's Sr Resin, Lot # SRS01127 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	136
Gravimetric Recovery	$\leq 105\%$	104 %
Chemical Recovery, Corrected	90-110 %	94.7 %



Jill Bryant
Quality Control Chemist

Prepared: January 23, 2007