

Certificate of Analysis

No. CA 13049

Product: Eichrom's Sr Resin


Lot #: SRS030813

Eichrom's Sr Resin Lot # SRS030813 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, Sr in 3.1 N nitric acid	≥ 100	138.5
% Gravimetric Recovery	≤ 105	99.4
% Chemical Recovery, Corrected	90-110	97.3
% Sr stripped in 0.05M HNO ₃		
Least acceptable mean (5 rep)	85.0	93.74
Max. standard deviation	2.5	1.4
%Ca stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.15	0.00
Max. standard deviation	0.05	0.00
%Ba stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.17
Max. standard deviation	0.10	0.07
%Y stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.15	0.0
Max. standard deviation	0.05	0.0
%Y in 8M HNO ₃ load and rinse		
Least acceptable mean (5 rep)	99.0	100.0
Max. standard deviation	2.5	0.0

Jill Bryant
Quality System Coordinator



April 11, 2013