

Certificate of Analysis

No. CA06042


Product: Eichrom's Sr Resin

Lot #: SRS03156

Eichrom's Sr Resin, Lot # SRS03156 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	≥ 100	122
Gravimetric Recovery	$\leq 105\%$	86.4 %
Chemical Recovery, Corrected	90-110 %	93.3 %



Jill Bryant

Quality Control Chemist

Prepared: April 14, 2006