



Certificate of Analysis

No. CA 15113

Product: Eichrom's Sr Resin

Lot #: SRS061115

Eichrom's Sr Resin Lot # SRS061115 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, Sr in 3.1 N nitric acid	≥ 100	121.5
% Gravimetric Recovery	≤ 105	100.2
% Chemical Recovery, Corrected	90-110	96.8
% Sr stripped in 0.05M HNO ₃		
Least acceptable mean (5 rep)	85.0	93.3
Max. standard deviation	2.5	1.7
%Ca stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.15	0.03
Max. standard deviation	0.05	0.01
%Ba stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.25	0.2
Max. standard deviation	0.10	0.1
%Y stripped in 0.05M HNO ₃		
Maximum acceptable mean (5 rep)	0.15	0.0
Max. standard deviation	0.05	0.0
%Y in 8M HNO ₃ load and rinse		
Least acceptable mean (5 rep)	99.0	100.0
Max. standard deviation	2.5	0.0

Sarah Teichma
Quality Manager

17 Jun 15

June 17, 2015