

# Certificate of Analysis

No. CA07152

Product: Eichrom's Sr Resin  
Lot #: SRS08107

Eichrom's Sr Resin Lot # SRS08107 has been manufactured according to Eichrom standard operating procedures. The resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

## ANALYSIS RESULTS

| <u>Test Value</u>                     | <u>Specification</u> | <u>Result</u> |
|---------------------------------------|----------------------|---------------|
| Appearance                            | White beads          | Pass          |
| Distribution Ratio, 3.1 N nitric acid | $\geq 100$           | 116           |
| Gravimetric Recovery                  | $\leq 105\%$         | 102 %         |
| Chemical Recovery, Corrected          | 90-110 %             | 95.5 %        |

  
\_\_\_\_\_  
Jill Bryant  
Quality Control Chemist

Prepared: October 30, 2007