

## Certificate of Analysis

No. CA05128

Product: Eichrom's Sr Resin


Lot #: SRS12055

Eichrom's Sr Resin, Lot # SRS12055 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	125
Gravimetric Recovery	$\leq 105\%$	100 %
Chemical Recovery, Corrected	90-110 %	95.9 %

Jill Bernstein



Quality Control Chemist

Prepared: December 20, 2005