

## Certificate of Analysis


No. CA06162

Product: Eichrom's Sr Resin Cartridges  
Lot #: SRSR6A

Eichrom's Sr Resin Cartridges, Lot # SRSR6A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets or exceeds all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	$\geq 100$	125
Gravimetric Recovery	$\leq 105\%$	100 %
Chemical Recovery, Corrected	90-110 %	95.9 %



Jill Bryant  
Quality Control Chemist

Prepared: December 20, 2006